Instruction for Use

Mucus Extractors

Pennine's Mucus Extractor range of devices are single use, covered plastic receptacle and associated tubing, containing no additives, intended to be used for the aspiration, and/or collection of mucus/sputum specimens required for further examination. To be used on all patient populations, by clinician or suitably trained person, in a clinical setting.

Instructions for Use:

Caution: FOR SINGLE USE, DO NOT REUSE OR RE-STERILISE, STERILE DEVICE.

Mucus Extractor with integral cannula (Codes pre-fixed ME-)
1. Select device with cannula sizing to suit the patient. 2. Remove Mucus extractor from packaging. 3. Check all contents are present and lid securely attached to the collection bottle. 4. Connect the mucus extractor to an appropriate length of suction connecting tubing. 5. With the suction source turned off, securely connect the suction tube to the suction source, turning on prior to inserting the cannula into the nasal oral cavity. 6. If using device with vacuum control, use fingertip to cover vent to control suction. 7. Follow local policies and procedures for mucus extraction and specimen collection. 8. Remove the suction cannula immediately upon collection of the sample. 9. When procedure is complete, remove the screw top lid/tubing, and replace with the spare screw cap. Ensure lid is secure. 10. The specimen can then be used for further examination or disposal.

Mucus Extractor without integral cannula (Codes pre-fixed MST-)
1. Remove Mucus extractor from packaging. 2. Check all contents are present and lid securely attached to the collection bottle. 3. Push fit the chosen suction cannula/catheter or device into the female connector. 4. Connect the mucus extractor to an appropriate length of suction connecting tubing. 5. With the suction source turned off, securely connect the suction tube to the suction source, turning on prior to inserting the cannula into the nasal oral cavity. 6. Follow local policies and procedures for mucus extraction and specimen collection. 7. Remove the suction cannula immediately upon collection of the sample. 8. When procedure is complete, remove the screw top lid/tubing, and replace with the spare screw cap. Ensure lid is secure. 9. The specimen can then be used for further examination or disposal.

Warnings:

EN

- Do not insert the catheter past the back of the nose or throat.
- Do not use if packaging is damaged or open, sterility will be compromised, there may be a risk of infection and the device may not perform as intended.
- Do not use if you observe any damage to, blockage of, or foreign bodies in the device as the device may not perform as intended or cause injury to the patient, user or other.
- Do not use if you observe any kinking in the tubing as this may reduce suction and the device may not perform as intended.
- Graduations on specimen bottles are for indication purposes only. The device does not have a measuring function.
- Device is intended for specimen collection only. Once sample is obtained remove the suction catheter/cannula from the patient.
- Device is not intended to preserve samples. Do not use the device to store samples.
- Do not reuse device as there may be a risk of infection to patient, user or other and the device may not perform as intended.
- Device is intended for single patient use.
- Follow local policies and procedures to store and transport samples for further examination.
- Do not use the device past its expiry date. As it may impact the performance of the device.
- If you experience problems using the device, contact your healthcare professional for advice. Serious incidents that occur in relation to the devices intended use, should be reported to the manufacturer and local competent authority.

Storage and Disposal:

- Store in a clean, cool, dry place away from direct sunlight.
- Dispose of as clinical waste in-line with local and Government guidelines.

DEHP FREE

Ivor Shaw Ltd Trading as Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY
Tel: +44 (0)1332 794880 Fax: +44 (0)1332 794890
www.penninehealthcare.co.uk

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