



# ENFit Nasogastric and Feeding Tube

With a unique enteral-specific design

# ENFit Nasogastric and Feeding Tube

The ENFit connector has a unique enteral-specific design that:

Does not allow connectivity with connectors for any other clinical use.

- Provides a locking feature for a secure, leak free connection
- Administration sets and syringes have (a) female connector ends that fit onto or over (a) male patient access feeding tube ports

The purpose of the new connector is to help reduce the risk of enteral tube feeding misconnections and improve patient safety. The new ISO standard, ISO-80369- 3, has been established for Luer connectors on the nutrition formula end and the patient-access end. There is just one standard connector that will be utilized by all feeding set manufacturers and universally adopted into practice.

Designed to promote improved patient safety and help ensure that connectors for unrelated delivery systems are incompatible.

Pennine will have the ENFit connector available to our customers on its Infant Feeding Tubes, Ryles Tube and Levins Tube product ranges.

**Please contact Pennine to order or for further information:**

Tel: +44 (0)1332 794880

[marketing@penninehealthcare.co.uk](mailto:marketing@penninehealthcare.co.uk)

[www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)

## Ryles / Feeding Nasogastric Tube

X-Ray Graduated with ENFit Connector, metal tip and four side eyes.

Code	Description	Box Qty
RT-2108/EF	8fg x 122cm	100 (10 x 10)
RT-2110/EF	10fg x 122cm	100 (10 x 10)
RT-2112/EF	12fg x 122cm	100 (10 x 10)
RT-2114/EF	14fg x 122cm	100 (10 x 10)
RT-2116/EF	16fg x 122cm	100 (10 x 10)
RT-2118/EF	18fg x 122cm	100 (10 x 10)

## Ryles / Feeding Nasogastric Tube

X-Ray Graduated with ENFit Connector, opaque closed tip and four side eyes.

Code	Description	Box Qty
RT-2306/105/EF	6fg x 105cm	100
RT-2306/EF	6fg x 122cm	100
RT-2308/EF	8fg x 122cm	100
RT-2310/EF	10fg x 122cm	100
RT-2312/EF	12fg x 122cm	100
RT-2314/EF	14fg x 122cm	100
RT-2316/EF	16fg x 122cm	100
RT-2318/EF	18fg x 122cm	100
RT-2320/EF	20fg x 122cm	50

## Levins / Feeding Nasogastric Tube

X-Ray Graduated with ENFit Connector, open tip and four side eyes.

Code	Description	Box Qty
XLT-2205/105/EF	5fg x 105cm	100
XLT-2206/EF	6fg x 122cm	100
XLT-2208/EF	8fg x 122cm	100
XLT-2210/EF	10fg x 122cm	100
XLT-2212/EF	12fg x 122cm	100
XLT-2214/EF	14fg x 122cm	100
XLT-2216/EF	16fg x 122cm	100
XLT-2218/EF	18fg x 122cm	100
XLT-2220/EF	20fg x 122cm	50

## Infant / Feeding Nasogastric Tube

X-Ray Graduated with ENFit Connector, closed tip and two side eyes.

Code	Description	Box Qty
FT-1905/40/EF	5fg x 40cm	100
FT-1906/40/EF	6fg x 40cm	100
FT-1908/40/EF	8fg x 40cm	100
FT-1910/40/EF	10fg x 40cm	100
FT-1905/60/EF	5fg x 60cm	100
FT-1906/60/EF	6fg x 60cm	100
FT-1908/60/EF	8fg x 60cm	100
FT-1910/60/EF	10fg x 60cm	100
FT-1905/80/EF	5fg x 80cm	100
FT-1906/80/EF	6fg x 80cm	100
FT-1908/80/EF	8fg x 80cm	100
FT-1910/80/EF	10fg x 80cm	100