

# Urology

Ur

Catheters

Pre-lubricated Catheters

Self-Catheterisation System

Urine Collection

# 02.

---

<b>Contents</b>		<b>Page</b>
<b>Catheters</b>	Nelaton Male Catheter	3
	Female Catheter	3
	Tiemann Catheter	4
<hr/>		
<b>Urine Collection</b>	Mid-Stream Urine Specimen Collector	4
<hr/>		
<b>Pre-lubricated Catheters</b>	PreLube (male)	5
	PreLube (female)	5
	PreLube Bag for PreLube Catheter (male)	6
	PreLube Tieman Catheter (male)	6
<hr/>		
<b>Self Catheterisation Systems</b>	Pennine Intro-Uri-Cath	7



## Nelaton (Male) Catheters

For intermittent bladder drainage.

### Nelaton (Male) Catheters

- Medical grade PVC
- Sterile, individual peel pouches
- 100 units per inner polythene bag
- 38cm length
- Soft rounded tip for atraumatic insertion
- Two side eyes for maximum drainage
- Supplied single wrapped (SW) as standard

Single Wrapped Product Code	Qty Per Box	Drug Tariff Code	Qty Per Box	Double Wrapped Product Code	Qty Per Box	CH/Diam
NC-1206/SW	500			NC-1206/DW	400	6ch ø 2.00mm
NC-1208/SW	500			NC-1208/DW	400	8ch ø 2.66mm
NC-1210/SW	500			NC-1210/DW	400	10ch ø 3.33mm
NC-1212/SW	500	NC-1212/FP	100	NC-1212/DW	400	12ch ø 4.66mm
NC-1214/SW	500	NC-1214/FP	200	NC-1214/DW	400	14ch ø 4.66mm
NC-1216/SW	500	NC-1216/FP	200	NC-1216/DW	400	16ch ø 5.33mm
NC-1218/SW	500			NC-1218/DW	400	18ch ø 6.00mm
NC-1220/SW	400			NC-1220/DW	300	20ch ø 6.66mm
NC-1222/SW	400			NC-1222/DW	300	22ch ø 7.33mm
NC-1224/SW	400			NC-1224/DW	300	24ch ø 8.00mm

### NHS supply chain codes

- Double wrapped (DW) to provide aseptic non-touch technique
- Ward pack configuration
- Pouched in inners of 25

Product Code	NHS NPC Code	Qty Per Box	CH/Diam
NC-1210/DW/TRE	FST627	400	10ch ø 3.33mm
NC-1212/DW/TRE	FST628	400	12ch ø 4.00mm
NC-1214/DW/TRE	FST629	400	14ch ø 4.66mm



## Female Catheters

For intermittent bladder drainage.

### Female Catheters

- 23cm length
- Medical grade PVC
- Sterile, individual peel pouches
- 100 units per inner polythene bag
- Soft rounded tip for atraumatic insertion
- Two side eyes for maximum drainage

Single Wrapped Product Code	Qty Per Box	Drug Tariff Code	Qty Per Box	CH/Diam
FC-1406	1000			6ch ø 2.00mm
FC-1408	1000			8ch ø 2.66mm
FC-1410	1000	FC-1410/FP	50	10ch ø 3.33mm
FC-1412	1000	FC-1412/FP	50	12ch ø 4.00mm
FC-1414	1000	FC-1414/FP	50	14ch ø 4.66mm
FC-1416	1000			16ch ø 5.33mm

### NHS supply chain codes

- Ward pack configuration
- Pouched in inners of 25
- 16 bags contain 25 units

Product Code	NHS NPC Code	Qty Per Box	CH/Diam
FC-1406/TRE	FST626	400	6ch ø 2.00mm
FC-1408/TRE	FST022	400	8ch ø 2.66mm
FC-1410/TRE	FST023	400	10ch ø 3.33mm
FC-1412/TRE	FST024	400	12ch ø 4.00mm

# 04.



## Tiemann Catheter

Traditional Nelaton with angled, pointed tip for directional stability.

- Slightly hardened tip for ease of insertion through areas of obstruction
- Medical grade PVC
- Sterile, individual peel packs, single wrapped
- Single use only
- 38cm length

Product Code	Qty Per Box	CH/Diam
NC-1708	100	8ch ø 2.66mm
NC-1710	100	10ch ø 3.33mm
NC-1712	100	12ch ø 4.00mm
NC-1714	100	14ch ø 4.66mm
NC-1716	100	16ch ø 5.33mm
NC-1718	100	18ch ø 6.00mm



## Mid-Stream Urine Specimen Collector

For aseptic collection of midstream urine.

- Easy to use, 30ml container
- Patient details label
- Patient instructions

Product Code	Type	Trap Size	Qty Per Box
MSU-1000	Without cleansing wipe	30ml	50
MSU-1001	With cleansing wipe	30ml	50

## PreLube Catheters

The Pennine PreLube is a pre-lubricated catheter for Intermittent Urethral Self-Catheterisation. The PreLube is coated with Lubrajel® to reduce friction and minimise discomfort for the patient. The unique coating eradicates the need for water or any other gelling agent. The catheter is completely enclosed, making the PreLube very safe to use because there is no need to touch the catheter itself. The PreLube is ideal for teaching self-catheterisation.



### PreLube (male)

Male PreLube (43cm)

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel packs, double wrapped
- Sealed pre-cut polythene wrapper
- Single use only
- Smooth closed end, two staggered eyes
- Hard PVC capsule containing Lubrajel® to offer total pre-lubrication of the Nelaton catheter and ensure quick and comfortable insertion
- Detachable packaging tip for aseptic introduction of the catheter into the urethra

Product Code	Qty Per Box	Drug Tariff Code	CH/Diam
NC-1606	100	–	6ch ø 2.00mm
NC-1608	100	NC-1608/FP	8ch ø 2.66mm
NC-1610	100	NC-1610/FP	10ch ø 3.33mm
NC-1612	100	NC-1612/FP	12ch ø 4.00mm
NC-1614	100	NC-1614/FP	14ch ø 4.66mm
NC-1616	100	NC-1616/FP	16ch ø 5.33mm



### PreLube (female)

Female PreLube (23cm)

Product Code	Qty Per Box	Drug Tariff Code	Qty Per Box	CH/Diam
FC-1606	100	–	–	6ch ø 2.00mm
FC-1608	100	FC-1608/FP	10	8ch ø 2.66mm
FC-1610	100	FC-1610/FP	10	10ch ø 3.33mm
FC-1612	100	FC-1612/FP	10	12ch ø 4.00mm
FC-1614	100	FC-1614/FP	10	14ch ø 4.66mm
FC-1616	100	FC-1616/FP	10	16ch ø 5.33mm

PreLube and PreLube Bags are available in the U.K. on Drug Tariff.

# 06.



## PreLube Bag for PreLube Catheter (male)

- 1 litre collection bag
- 43cm Nelaton catheter
- Easily disposable

Product Code	Qty Per Box	Drug Tariff Code	CH/Diam
NC-1608/BG	100	NC-1608/BG/FP	8ch ø 2.66mm
NC-1610/BG	100	NC-1610/BG/FP	10ch ø 3.33mm
NC-1612/BG	100	NC-1612/BG/FP	12ch ø 4.00mm
NC-1614/BG	100	NC-1614/BG/FP	14ch ø 4.66mm
NC-1616/BG	100	NC-1616/BG/FP	16ch ø 5.33mm

## PreLube Tieman Catheter (male)

PreLube catheter with angled pointed tip for directional stability

Product Code	Qty Per Box	CH/Diam
NC-1908	100	8ch ø 2.66mm
NC-1910	100	10ch ø 3.33mm
NC-1912	100	12ch ø 4.00mm
NC-1914	100	14ch ø 4.66mm
NC-1916	100	16ch ø 5.33mm

PreLube and PreLube Bags are available in the U.K. on Drug Tariff.



## Pennine Intro-Uri-Cath

Sterile Field Urinary Catheter Introducer. A complete intermittent urethral catheterisation procedure kit. Extends the sterile field into the patient.

- Easy, quick, time saving, economical, safe evidence of reduced bacteriuria
- Sterile, individual peel pack with clear instructions for use
- Standard Nelaton catheter
- Silicone introducer with lubricant
- Graduated bag for dual purpose no-touch technique and collection of urine

Product Code	Qty Per Box	CH/Diam
NC-1508	20	8ch ø 2.66mm
NC-1510	20	10ch ø 3.33mm
NC-1512	20	12ch ø 4.00mm
NC-1514	20	14ch ø 4.66mm



## Indications

### Hospitals

- Where intermittent catheterisation is preferable to use of indwelling catheters
- Post trauma and post surgical bladder relief
- Prior to lower abdominal surgery
- Acute retention of urine and all situations where there is inadequate bladder drainage
- Spinal injuries, Spina Bifida where there is risk of U.T.I.
- Teaching self catheterisation
- Elderly and other patient groups where there is a susceptibility to nosocomial U.T.I.

### Maternity

- Pre and post delivery bladder drainage to reduce chances of bacteriuria
- Measurement of residual urine

### Rehabilitation/Community

- For teaching self catheterisation
- For many patient groups where long term intermittent catheterisation is associated with a need for a complete, self contained system and which reduces the risk of induced U.T.I.
- For general practitioner and other community medical personnel for emergency use, e.g. urinary retention

**Pennine**

City Gate, London Road, Derby, DE24 8WY

Tel: +44 (0)1332 794880

[marketing@penninehealthcare.co.uk](mailto:marketing@penninehealthcare.co.uk)

[www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)