

CONSUMABLES



01.

Laparoscopic Surgery

A core range of product specifically designed for use in laparoscopic surgical procedures.

Veress Needle

Intended for the safe, simple and effective creation of the pneumoperitoneum.

Product Code	Description
AELAP0001	120mm Veress Needle
AELAP0010	150mm Veress Needle



Camera Drape

The disposable camera sleeve is a sterile single-use drape made of tough polyethylene.

The drape is the quickest and simplest way to create a sterile field for use with video/camera systems during operative laparoscopic procedure.

Product Code	Description
E6222	Camera cover



Consumables

Laparoscopic Surgery

Insufflation Tubing Sets

For use in keyhole or minimally invasive surgery. Intended as a connection and filtration device between a pressure regulator carbon dioxide insufflator and a Veress needle. Designed to reduce the risks of backflow of body fluids into insufflator.

Product Code	Length (m)	Inside Dia (mm)	Product Description	Qty Per Box
AEINS1250	2.5	6	Hydrophobic filter, Male Luer/Wolf Connector	10
AEINS1330E	3.0	6	Hydrophobic filter, Open end/Funnel Connector	10
AEINS1330LL	3.0	6	Hydrophobic filter, Male Luer/Male Luer	10
AEINS1331	3.0	6	Hydrophobic filter, Male Luer/15/22mm Filter Adaptor	10
AEINS1333	3.0	6	Hydrophobic filter, Male Luer/Funnel Connector	10
AEINS1334	3.0	6	Hydrophobic filter, Female Luer/Male Luer	10

AEINS1333



AEINS1330LL



AEINS1330E



AEINS1250



AEINS1334



AEINS1331



Esophageal Stethoscope

Used to ascertain the condition of the heart and lungs by listening to their action.

We have short (55cm) and long (110cm) sizes available.

55cm Esophageal Stethoscope

Product Code	Description
ESV-5024	55cm Esophageal Stethoscope fg24
ESV-5018	55cm Esophageal Stethoscope fg18
ESV-5012	55cm Esophageal Stethoscope fg12



110cm Esophageal Stethoscope

Product Code	Description
ESV-1124	110cm Esophageal Stethoscope fg24
ESV-1118	110cm Esophageal Stethoscope fg18
ESV-1112	110cm Esophageal Stethoscope fg12

Frazier Suction Probes

Single use frazier probe manufactured from anodized aluminium for extra flexibility. With 'VC' vacuum control hole for precise control over the amount of suctioning. Used in open surgical procedures where delicate suctioning is required.

Product Code	Length (cm)	Size (ch)
AEMF1108	23	8
AEMF1110	23	10
AEMF1112	23	12



Microflex Suction Probes

PVC disposable suction probes suitable for intricate, controlled suction.


Product Code	Length (cm)	Size (mm)
YS-3088	15	2.6
YS-3108	15	3.33
YS-3128	15	4.0



Feeding Tubes

48cm tube for enteral feeding. Closed end, two side eyes. Supplied in sterile individual peel pouches.

Product Code	Length (cm)	Size (mm)
FT-1605/48	48	1.66
FT-1606/48	48	2.00
FT-1608/48	48	2.66
FT-1610/48	48	3.33



Chemotherapy Products

Chemo Apparel

Product Code	Description
AIG4016	Extra-large Gown Isolation/Procedure
VT0034	Chemo Face Mask

Chemo Prep Pack

Product Code	Description
VKT-08800	Chemo Prep Pack
1	Blue wrap 60cm x 60cm
1	Gloves, powder free, medium, Walleted
1	Chemo protection pad
1	Chemo spike, with venting filter, 0.2um
1	Yellow Waste bag
2	Wipe, 70% alcohol, 2% chlorhexidine gluconate, PDI
1	Face mask with visor
2	Oversleeves with elastic

Surgical Consumables

Cautery / Diathermy Pencil and Tip Cleaner

A single use, disposable, mono-polar electrosurgical pencil for soft tissue surgery. 3m connecting cable and a 2.5cm removable blade. Tip cleaner; for the removal of cauterised tissue from the pencil. Sterile, disposable and adhesive backed for attachment to drape or trolley cover. Supplied in sterile peel pouches.



Product Code	Description
SW12200	Diathermy Pencil
SW12200-B10	Diathermy Pencil set with tip cleaner and holster
TP-4711/S	Disposable Adson Bipolar forceps with flying lead 15mm and 0.5 tip

06.

Apparel

Product Code	Description
AFM2003	Surgical Masks

A range of gloves available upon request

Pennine 



Pennine Veterinary
Ivor Shaw Limited t/a Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, UK.
Tel: +44 (0)1332 794877
www.pennineveterinary.co.uk