

# AIRWAY MANAGEMENT

---





# CONTENTS

About Pennine	3
Pennine Airway Management Range	4

Catheters	
Funnel Suction Catheters	5
Oppo-Cath Suction Catheters	6
Prestrol Suction Catheters	7
Prestrol Oppo-Cath Suction Catheters	8
Clear Graduated Oppo-Cath Suction Catheters	9
Angled Tip Suction Catheters	10

Mucus Extractors	
Mucus Specimen Traps	11
Mucus Extractors	12

Oxygen Therapy	
Guedel Airways	13
Link Oxygen Connecting Tubing	14



## ABOUT PENNINE

Pennine Healthcare is proud to be one of the UK's leading, employee owned, manufacturer and distributor, of single use, sterile Ward & Theatre Consumables and Custom Procedure Packs, and since our incorporation in 1963, the business has grown year on year. We market our products and services in over 50 countries worldwide and employ around 200 people at our purpose-built premises in Derby.

In line with our ISO 13485 accreditation, Pennine has a policy of continual improvement in order to meet the needs of the changing global standards for sterile, single use medical devices. Continual investment in people, training, and technology, has enabled us to provide an increasingly wide range of user friendly, high specification products. We also promote standards of safety, health, and welfare to provide a healthy and safe working environment.

Our vision will be delivered by supplying world class products and delivering distinctive customer service supported by a culture of continuous

improvement in investing in our people. Our core skills in manufacturing and engineering make us competitive through technological innovation, enabling us to be innovative, flexible and responsive, supplying healthcare organisations and businesses across the globe with the highest quality British manufactured products.

Pennine Healthcare is also proud to be one of the leading suppliers of PPE into healthcare settings, and as a hospital focussed company, we sell directly to healthcare professionals in all areas, including intensive care and operating theatres. In 2014 Pennine Healthcare saw the need to develop a range of PPE and Infection Control products, due to a lack of complete, high quality products available for the general public, and thus ARMA was born.

We hope you find what you are looking for in the following pages of this brochure, but if you do not, please do not hesitate to get in touch. We're a friendly bunch and love speaking to new, and old, customers alike.



# PENNINE AIRWAY MANAGEMENT RANGE

Our Airways Management Range offers superb quality, performance, and solution, to your airway management needs.

We strive to achieve the following in our Airways Management Range:

## CHOICE


Whether you are looking for Catheters, Guedel Airways, or more, our product range involves a variety of lengths, diameters, and connectors, to fulfil yours and your patients needs and requirements.



## SPEED


Our Airways Management Range are colour coded, for easy size identification. In critical emergencies the colour coding\* system allows for instant identification of the size of the product contained within the pouch.

\* Colour coding of products comply to ISO 5364:2016






## SAFETY

All of our products are Latex and DEHP free as standard, putting patient safety and the forefront of everything we supply.



In line with our ISO 13485 accreditation, Pennine has a policy of continual improvement in order to meet the needs of the changing international standards for sterile, single use medical devices.

# FUNNEL SUCTION CATHETERS



## 48cm/53cm/60cm

The product is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient. Critically ill or injured patients who require a breathing (endotracheal) tube or tracheostomy tube need occasional suctioning to remove secretions from the airway.

- Medical grade PVC
- Sterile individual peel pouches
- Colour coded funnel end
- Length excludes connector
- Standard configuration= terminal eye, 1 side eye
- Double offset eye option=/DE
- Terminal eye only option=/TO
- 100 per inner box

48cm Single Eye	48cm Double Eye	48cm Terminal Eye Only	Qty Per Box	CH/Diam	
* AESC100648	AESC100648DE	AESC100648TO	100	6ch ø 2.0mm	
* AESC100848	AESC100848DE	AESC100848TO	100	8ch ø 2.7mm	
* AESC101048	AESC101048DE	AESC101048TO	100	10ch ø 3.3mm	
* AESC101248	AESC101248DE	AESC101248TO	100	12ch ø 4.0mm	
* AESC101448	AESC101448DE	AESC101448TO	100	14ch ø 4.7mm	
* AESC101648	AESC101648DE	AESC101648TO	100	16ch ø 5.3mm	
* AESC101848	AESC101848DE	AESC101848TO	100	18ch ø 6.0mm	
* AESC102048	-	-	100	20ch ø 6.7mm	


53cm Single Eye	53cm Double Eye	60cm Single Eye	Qty Per Box	CH/Diam	
* AESC100653	AESC100653DE	-	100	6ch ø 2.0mm	
* AESC100853	AESC100853DE	-	100	8ch ø 2.7mm	
* AESC101053	AESC101053DE	AESC101060	100	10ch ø 3.3mm	
* AESC101253	AESC101253DE	AESC101260	100	12ch ø 4.0mm	
* AESC101453	AESC101453DE	AESC101460	100	14ch ø 4.7mm	
* AESC101653	AESC101653DE	AESC101660	100	16ch ø 5.3mm	
* AESC101853	AESC101853DE	AESC101860	100	18ch ø 6.0mm	



Want more information?  
Get in touch on +44 (0) 1332 794880 or email [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)

Please note: All dimensions are nominal.  
Please note any products featuring a "\*" indicates the product is distributed by Pennine.  
For Regulatory and Technical information, Technical Data Sheets, are available upon request.

OPPO-CATH SUCTION CATHETERS




53cm

The product is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient. Critically ill or injured patients who require a breathing (endotracheal) tube or tracheostomy tube need occasional suctioning to remove secretions from the airway.

- Colour coded funnel end
- Sterile individual peel pouches
- Terminal eye, two small opposed side eyes
- 100 per inner box

53cm Length	Qty Per Box	CH/Diam	
* AEOSC130653	100	6ch ø 2.0mm	
* AEOSC130853	100	8ch ø 2.7mm	
* AEOSC131053	100	10ch ø 3.3mm	
* AEOSC131253	100	12ch ø 4.0mm	
* AEOSC131453	100	14ch ø 4.7mm	
* AEOSC131653	100	16ch ø 5.3mm	
* AEOSC131853	100	18ch ø 6.0mm	

PRESTROL SUCTION CATHETERS



48cm/60cm

The product is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient. Critically ill or injured patients who require a breathing (endotracheal) tube or tracheostomy tube need occasional suctioning to remove secretions from the airway.

- Medical grade PVC
- Sterile individual peel pouches
- Colour coded Prestrol end enabling vacuum control
- Standard configuration= terminal eye, 1 side eye
- 400 per inner box

48cm Length	60cm Length	Qty Per Box	CH/Diam	
* SC-1105/48	SC-1105/60	400	5ch ø 1.7mm	
* SC-1106/48	SC-1106/60	400	6ch ø 2.0mm	
* SC-1108/48	SC-1108/60	400	8ch ø 2.7mm	
* SC-1110/48	SC-1110/60	400	10ch ø 3.3mm	
* SC-1112/48	SC-1112/60	400	12ch ø 4.0mm	
* SC-1114/48	SC-1114/60	400	14ch ø 4.7mm	
* SC-1116/48	SC-1116/60	400	16ch ø 5.3mm	
* SC-1118/48	SC-1118/60	400	18ch ø 6.0mm	



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)

PRESTROL OPPO-CATH SUCTION CATHETERS




48cm/53cm/60cm

The product is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient. Critically ill or injured patients who require a breathing (endotracheal) tube or tracheostomy tube need occasional suctioning to remove secretions from the airway.

- Medical grade PVC
- Sterile individual peel pouches
- With colour coded Prestrol Vacuum Control Device
- Terminal eye, two small opposed side eyes
- 400 per inner box

48cm Length	53cm Length	60cm Length	Qty Per Box	CH/Diam	
* OSC-1505/48	OSC-1505/53	OSC-1505/60	400	5ch ø 1.7mm	
* OSC-1506/48	OSC-1506/53	OSC-1506/60	400	6ch ø 2.0mm	
* OSC-1508/48	OSC-1508/53	OSC-1508/60	400	8ch ø 2.7mm	
* OSC-1510/48	OSC-1510/53	OSC-1510/60	400	10ch ø 3.3mm	
* OSC-1512/48	OSC-1512/53	OSC-1512/60	400	12ch ø 4.0mm	
* OSC-1514/48	OSC-1514/53	OSC-1514/60	400	14ch ø 4.7mm	
* OSC-1516/48	OSC-1516/53	OSC-1516/60	400	16ch ø 5.3mm	
* OSC-1518/48	OSC-1518/53	OSC-1518/60	400	18ch ø 6.0mm	

CLEAR GRADUATED OPPO-CATH SUCTION CATHETERS



35cm/48cm/60cm

The product is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient. Critically ill or injured patients who require a breathing (endotracheal) tube or tracheostomy tube need occasional suctioning to remove secretions from the airway.

- Medical grade PVC
- Clear Prestrol Vacuum control
- Graduations every 1cm
- Terminal eye, two small opposed side eyes
- Sterile individual peel pouches
- 400 per inner box

35cm Length	48cm Length	60cm Length	Qty Per Box	CH/Diam	
* OSC-1205/35	OSC-1205/48	OSC-1205/60	400	5ch ø 1.7mm	
* OSC-1206/35	OSC-1206/48	OSC-1206/60	400	6ch ø 2.0mm	
* OSC-1207/35	OSC-1207/48	OSC-1207/60	400	7cm ø 3.0mm	
* OSC-1208/35	OSC-1208/48	OSC-1208/60	400	8ch ø 2.7mm	
* OSC-1210/35	OSC-1210/48	OSC-1210/60	400	10ch ø 3.3mm	
* OSC-1212/35	OSC-1212/48	OSC-1212/60	400	12ch ø 4.0mm	
* OSC-1214/35	OSC-1214/48	OSC-1214/60	400	14ch ø 4.7mm	
* OSC-1216/35	OSC-1216/48	OSC-1216/60	400	16ch ø 5.3mm	




Want more information?  
Get in touch on +44 (0) 1332 794880 or email [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)

ANGLED TIP SUCTION CATHETERS




53cm

The product is used to provide access to respiratory tract to remove secretion with maximum comfort to the patient. Critically ill or injured patients who require a breathing (endotracheal) tube or tracheostomy tube need occasional suctioning to remove secretions from the airway. The respiratory therapist, nurse, or ICU technician suctions a patient by inserting a small flexible tube called a suction catheter into the breathing tube.

- Extra smooth low friction surface tubing facilitates easy intubation
  - A traumatic, soft & rounded open tip with two lateral eyes
- Colour coded connector for easy size identification
  - Sterile, supplied in individual peel pouches

53cm Length		Qty Per Box	CH/Diam	
* AEOSC171053		100	10ch ø 3.3mm	
* AEOSC171253		100	12ch ø 4.0mm	
* AEOSC171453		100	14ch ø 4.7mm	
* AEOSC171653		100	16ch ø 5.3mm	

MUCUS SPECIMEN TRAPS



Our Mucus Specimen Traps are used to collect mucus specimens for further examinations and are useful in several situations.

- Sterile, individual pack
  - 25ml or 70ml trap
  - Attachment for suction catheter and suction connecting tube
- Patient details label
  - Spare cap, screw fit

Product Code	Capacity	Type	Qty Per Box
MST-3000	25ml	With spare screw cap - 14cm tubing	50
MST-3070	25ml	With reduced length connection to suction catheter to ease movement of secretions to the trap. With spare screw cap & patient label.	50
MST-3180	25ml	With spare screw cap - 29cm tubing	75
MST-4000	70ml	For use in bronchoscopy. - 11cm tubing. With spare screw cap	50
MST-5000	25ml	With spare screw cap - 9cm tubing with VC	50
MST-7000	70ml	With male connector	50




For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)




For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)

MUCUS EXTRACTORS

Standard Funnel Control Connector



Vacuum Control Connector



Our Mucus Extractors are devices used for aspiration of mucus or secretions from the upper airways.

- 20ml Trap
- With colour coded cap for easy identification of catheter size
- With standard funnel or vacuum control connector for attachment to suction source
- Spare cap, screw fit



Without Vacuum Control Connector

Product Code	Suction Tube Length	Qty Per Box	CH/Diam	
ME-8108	25cm	75	8ch ø 2.7mm	
ME-8110	25cm	75	10ch ø 3.3mm	
ME-8112	25cm	75	12ch ø 4.0mm	
ME-8114	25cm	75	14ch ø 4.7mm	

With Vacuum Control Connector

Product Code	Suction Tube Length	Qty Per Box	CH/Diam	
ME-9108	25cm	75	8ch ø 2.7mm	
ME-9110	25cm	75	10ch ø 3.3mm	
ME-9112	25cm	75	12ch ø 4.0mm	
ME-9114	25cm	75	14ch ø 4.7mm	

GUEDEL AIRWAYS



Airways is a medical device called an airway adjunct used to maintain a patient (open) airway. It prevents the tongue from covering the epiglottis, which could prevent the person from breathing.

- Sterile, individual peel pouches
- Manufactured in linear low-density polyethylene
- Colour-coded bite-blocks for easy size identification
- Complies to ISO 5364:2016

Product Code	Description	Length	Qty Per Box	ISO Size	
* GA-93000	Guedel Airway, 4cm	40mm	100	4	
* GA-9300	Guedel Airway, 5cm	50mm	100	5	
* GA-9310	Guedel Airway, 6cm	60mm	100	6	
* GA-9311	Guedel Airway, 7cm	70mm	100	7	
* GA-9312	Guedel Airway, 9cm	90mm	80	9	
* GA-9313	Guedel Airway, 10cm	100mm	100	10	
* GA-9314	Guedel Airway, 11cm	110mm	100	11	
* GA-9315	Guedel Airway, 12cm	120mm	100	12	





Want more information?  
Get in touch on +44 (0) 1332 794880 or email [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)



# OXYGEN ADMINISTRATION TUBING



Pennine Oxygen administration tubes are non-sterile, single use devices intended to provide an additional length of tubing between the oxygen source and patient oxygen delivery device, in patients requiring supplementary oxygen from adults to paediatrics in both clinical, ambulatory, and home care environments. To be used by trained clinicians / homecare user, in clinical, homecare or emergency settings.

- Clean, non-sterile
- Individual polythene packs
- Also available in 30m rolls
- Link female, funnel ends can be easily cut to provide sure fit to any oxygen therapy device

### Standard Tubing

Product Code	Tube Length	Qty Per Box
OC-7311	1.8m	100
OC-7341	2.7m	100

### Rolls in Dispenser Boxes

Product Code	Tube Length	Qty Per Box
OL-8400	30m Coil	15
OL-8450	50m Coil	12



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to [customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)

Pennine



Ivor Shaw Limited t/a Pennine Healthcare  
300 City Gate, London Road, Derby, DE24 8WY, UK.  
Tel: +44 (0)1332 794880  
[www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)

@arma\_health Pennine Healthcare @Pennine\_Health