

NON-STERILE PRODUCT RANGE



CONTENTS

| | |
|--------------------------------------|----|
| About Pennine | 3 |
| Pennine Non-Sterile Range | 4 |
| Non-Sterile Link Yankauers | 5 |
| Non-Sterile Link Yankauer Sets | 6 |
| Non-Sterile Poole Suction Probes | 7 |
| Non-Sterile Poole Suction Probe Sets | 8 |
| Non-Sterile Micro-Flex Probes | 9 |
| Non-Sterile Suction Probes | 10 |
| Non-Sterile Suction Probe Sets | 11 |
| Non-Sterile Link Bubble Tubing | 12 |



ABOUT PENNINE

Pennine Healthcare is proud to be one of the UK's leading, employee owned, manufacturer and distributor, of single use, sterile and non-sterile Ward & Theatre Consumables and Custom Procedure Packs, and since our incorporation in 1963, the business has grown year on year. We market our products and services in over 50 countries worldwide and employ around 200 people at our purpose-built premises in Derby.

In line with our ISO 13485 accreditation, Pennine has a policy of continual improvement in order to meet the needs of the changing global standards for sterile, single use medical devices. Continual investment in people, training, and technology, has enabled us to provide an increasingly wide range of user friendly, high specification products. We also promote standards of safety, health, and welfare to provide a healthy and safe working environment.

Our vision will be delivered by supplying world class products and delivering distinctive customer service supported by a culture of continuous improvement in investing in our people. Our core

skills in manufacturing and engineering make us competitive through technological innovation, enabling us to be innovative, flexible and responsive, supplying healthcare organisations and businesses across the globe with the highest quality British manufactured products.

Pennine Healthcare is also proud to be one of the leading suppliers of PPE into healthcare settings, and as a hospital focussed company, we sell directly to healthcare professionals in all areas, including intensive care and operating theatres. In 2014 Pennine Healthcare saw the need to develop a range of PPE and Infection Control products, due to a lack of complete, high-quality products available for the general public, and thus ARMA was born.

We hope you find what you are looking for in the following pages of this brochure, but if you do not, please do not hesitate to get in touch. We're a friendly bunch and love speaking to new, and old, customers alike.



PENNINE NON-STERILE PRODUCT RANGE

Our Bulk Non-Sterile range offers superb quality, performance, and solution, to your surgical needs.

We strive to achieve the following in our Non-Sterile range:


CHOICE

Whether you are looking for Link Yankauers, Suction Probes, or more, our product range offers a variety of products in a range of lengths, diameters, and connectors, to fulfil yours and your patients' needs and requirements.




SPEED

Our Non-Sterile range provides quicker lead times, reduced packaging, and cheaper costs. Which means you will receive the devices you require, quicker, more cost effective and in a more sustainable manner.





SAFETY

We are fully MDR Certified. All of our products are Latex free and not manufactured using DEHP as standard, putting patient safety and the forefront of everything we manufacture.



We are delighted to announce that Pennine Healthcare have been awarded our MDR Certification for Class I Sterile, Class IIa and Procedure Packs.

NON-STERILE LINK YANKAUERS



Single use, rigid suction devices for use with suction connecting tubing for general purpose surgical suctioning.

- Available with or without vacuum control (VC) vent
- Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip



| Product Code | Type | Length | CH | Qty per Box |
|-----------------|--------------------------------|--------|------|-------------|
| LY-3400/BULK | Short Mini Tip | 190mm | 15ch | 300 |
| LY-3401/BULK | Short Mini Tip VC | 190mm | 15ch | 300 |
| LY-3500/BULK | Mini Fine Tip | 255mm | 12ch | 300 |
| LY-3501/BULK | Mini Fine Tip (VC) | 255mm | 12ch | 300 |
| LY-3600/BULK | Midi Tip | 240mm | 18ch | 300 |
| LY-3600/PR/BULK | Midi Tip (PR) | 240mm | 18ch | 300 |
| LY-3601/BULK | Midi Tip VC | 240mm | 18ch | 300 |
| LY-3700/BULK | Maxi Tip | 275mm | 22ch | 300 |
| LY-3701/BULK | Maxi Tip (VC) | 275mm | 22ch | 300 |
| LY-3800/BULK | Maxi Plus Tip | 280mm | 28ch | 300 |
| LY-3801/BULK | Maxi Plus Tip (VC) | 280mm | 28ch | 300 |
| LY-4400/BULK | Short Mini Super Tip | 195mm | 15ch | 300 |
| LY-4401/BULK | Short Mini Super Tip (VC) | 195mm | 15ch | 300 |
| LY-4420/BULK | Short Mini Super Tip (2 side) | 195mm | 15ch | 300 |
| LY-4500/BULK | Mini Fine Super Tip | 255mm | 12ch | 300 |
| LY-4500/NL/BULK | Mini Fine Super Tip (INT) | 255mm | 12ch | 300 |
| LY-4501/BULK | Mini Fine Super Tip (VC) | 255mm | 12ch | 300 |
| LY-4501/NL/BULK | Mini Fine Super Tip (VC) (INT) | 255mm | 12ch | 300 |
| LY-4600/BULK | Midi Super Tip | 240mm | 18ch | 300 |
| LY-4600/PR/BULK | Midi Super Tip (PR) | 240mm | 18ch | 300 |
| LY-4601/BULK | Midi Super Tip (VC) | 240mm | 18ch | 300 |
| LY-4700/BULK | Maxi Super Tip | 275mm | 22ch | 300 |
| LY-4700/PR/BULK | Maxi Super Tip (PR) | 275mm | 22ch | 300 |
| LY-4701/BULK | Maxi Super Tip (VC) | 275mm | 22ch | 300 |
| LY-3601/PC | Midi Tip VC | 240mm | 18ch | 50 |
| LY-3401/PC | Short Mini Tip VC | 190mm | 15ch | 50 |
| LY-3600/PC | Midi Tip | 240mm | 18ch | 50 |



Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk

Please note: All dimensions are nominal.
Please note any products featuring a "*" indicates the product is distributed by Pennine.
For Regulatory and Technical information, Technical Data Sheets, are available upon request.

NON-STERILE LINK YANKAUER SETS




Single use, Link Suction Yankauers, pre-connected to Pennine suction connecting tubing in 2m and 3metre lengths.

- Available with or without vacuum control (VC) vent
- Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip
- Tubing available in 2m and 3m lengths

| Product Code | Tubing Diam/Length | Qty per box | Yankauer Type |
|---------------|--------------------|-------------|--------------------------|
| LYS-5520/BULK | 6mm / 2m | 100 | LY-3500 Mini Fine Tip |
| LYS-5530/BULK | 6mm / 3m | 50 | LY-4400 Short Mini |
| LYS-5531/BULK | 6mm / 3m | 50 | LY-3501 Mini Fine Tip VC |
| LYS-5620/BULK | 6mm / 2m | 100 | LY-3600 Midi |
| LYS-5621/BULK | 6mm / 2m | 100 | LY-3601 Midi VC |
| LYS-5630/BULK | 6mm / 3m | 50 | LY-3600 Midi |
| LYS-5631/BULK | 6mm / 3m | 50 | LY-3601 Midi VC |
| LYS-5720/BULK | 6mm / 2m | 100 | LY-3700 Maxi |
| LYS-5721/BULK | 6mm / 2m | 100 | LY-3701 Maxi VC |
| LYS-5730/BULK | 6mm / 3m | 50 | LY-3700 Maxi |
| LYS-5731/BULK | 6mm / 3m | 50 | LY-3701 Maxi VC |
| LYS-6620/BULK | 7mm / 2.2m | 50 | LY-3600 Midi |
| LYS-6830/BULK | 7mm / 3m | 50 | LY-3800 Maxi Plus |
| LYS-7730/BULK | 8mm / 3m | 50 | LY-3700 Maxi |

NON-STERILE POOLE SUCTION PROBES



Single use, flexible sump suction probe with removable inner suction probe. Specifically for abdominal surgery, reducing trauma to surrounding tissue.

- Available with multiple side eyes to reduce tissue grab and for quick drainage and suctioning
- Removable Inner probe
- Available with straight or angled tip

| Product Code | Type | Length mm | CH/Tip | Qty per box |
|--------------|-----------------------|-----------|--------|-------------|
| YS-3006/BULK | Flexible straight tip | 300 | 20CH | 300 |
| YS-3008/BULK | Flexible (angled tip) | 300 | 20CH | 300 |





Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

NON-STERILE POOLE SUCTION PROBE SETS





Single use, disposable, suction probes pre attached to 2m and 3m connecting tube.

- Probes available 20ch and 30ch
- Pre-attached to 2m and 3m Tubing

| Product Code | Tubing Diam/Length | Qty per box | Poole Type |
|---------------|--------------------|-------------|---------------------------|
| SYS-5920/BULK | 6mm x 2m | 50 | YS-3006 Flexible Straight |
| SYS-5930/BULK | 6mm x 3m | 50 | YS-3006 Flexible Straight |

NON-STERILE MICRO-FLEX PROBES



Single use, suction probes with fingertip vacuum control. Designed for surgery and in applications where fine intricate controlled suction is required.

- Flexible suction tube with fingertip vacuum control
- Available in a range of lengths to suit a wide range of procedures.

| Product Code | Length | CH/Tip | Qty per box |
|-----------------|--------|---------------|-------------|
| YS-3128/BULK | 21.3cm | 12ch ø 4.00mm | 500 |
| YS-3148/14/BULK | 22.3cm | 14ch ø 4.66mm | 500 |




For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

NON-STERILE SUCTION PROBES




*SP-3800/BULK pictured

Single use, disposable, suction probes pre attached to 2m and 3m connecting tube.

- Probes available 20ch and 30ch
- Pre-attached to 2m and 3m Tubing

| Product Code | Type | Length | CH/Tip | Qty per Box |
|--------------|--|--------|--------|-------------|
| SP-3800/BULK | Suction probe with male connector | 245mm | 20CH | 300 |
| SP-3900/BULK | Flat Tip Suction Probe with male connector | 250mm | 20CH | 300 |
| SP-4200/BULK | Maxi Suction Probe | 280mm | 30CH | 300 |
| SP-4300/BULK | Flat Tip Maxi Suction Probe | 255mm | 30CH | 300 |
| YS-3000/BULK | Bent Tip Suction Probe | 265mm | 20CH | 300 |
| YS-3001/BULK | Bent Tip Suction Probe (VC) | 265mm | 20CH | 300 |

NON-STERILE SUCTION PROBE SETS



*SP-3800/BULK pictured

Single use, disposable, suction probes pre attached to 2m and 3m connecting tube.

- Probes available 20ch and 30ch
- Pre-attached to 2m and 3m Tubing

| Product Code | CH/Tip | Qty per Box | Tube ID | Probe Type |
|-----------------|--------|-------------|---------|-------------------|
| SPS-5020/W/BULK | 20CH | 100 | 6mm | Standard Probe |
| SPS-5021/BULK | 20CH | 100 | 6mm | Standard Probe VC |
| SPS-6020/BULK | 20CH | 100 | 8mm | Maxi Probe |





Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

NON-STERILE LINK BUBBLE TUBING




Link bubble tubing - available in 30 and 50 metre lengths supplied non-sterile in convenient dispenser bags.

- Tube extruded with a bubble at approx 1m intervals
- Simply cut through the bubble at the required length to create either male or female end which will fit a suction device
- Three diameters available: 5mm, 6mm and 7mm

| Product Code | Type | Length | CH/Tip | Inside Diam. | Qty per Box |
|-------------------------|-------|--------|--------|--------------|-------------|
| SL-8005 | Plain | 30m | 24CH | 5mm | 8 Coils |
| SL-8055 | Plain | 50m | 24CH | 5mm | 4 Coils |
| SL-8356 | Plain | 50m | 26CH | 6mm | 4 Coils |
| SL-8007 (International) | Plain | 30m | 28CH | 7mm | 7 Coils |

TUBING CONNECTOR RIGID




A short length rigid PVC tube that tapers at each end to form male connectors allowing two sections of suction connecting tubing to be joined to increase the length.

- Can be used to extend length of tubing (i.e to connect two tubes together of either the same or different diameters).

| Product Code | Qty per Box |
|--------------|-------------|
| DP-501/BULK | 500 |

CAPPED VACUUM CONTROL CONNECTOR



The vacuum control unit is comprised of two PVC male connectors with a hole in between that can be covered by a finger during use to control the vacuum and suction flow rate of a device while in use. The vacuum control unit also has a cap that can be placed to cover the hole to allow for continuous suction without the need of the user to cover the hole with their finger.

- Can be used to connect two devices and allow control of fluid flow. Male ends push fit into suction tubing and/or device. Open cap and use finger pressure to control the vacuum flow.

| Product Code | Qty per Box |
|---------------|-------------|
| VC-2/100/BULK | 100 |
| VC-2/BULK | 1200 |




For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

NON-STERILE SUCTION CONNECTING TUBE



Single Use, Suction connecting tubing to connect suction sources to suction waste collection systems.

- Secure fit to most suction devices and sources
- Available with optional male connector with or without VC

| Product Code | Type | Length | CH/Tip | Inside Diam. | Qty per Box |
|-----------------|----------------|--------|--------|--------------|-------------|
| CT-4012/NS | F/F | 1.8m | 26CH | 6mm | 25 |
| CT-4014/NS | F/F/M | 1.8m | 26CH | 6mm | 25 |
| CT-4020/BULK | Flush cut ends | 2m | 26CH | 6mm | 100 |
| CT-4022/BULK | F/F | 2m | 26CH | 6mm | 100 |
| CT-4024/BULK | F/F/M | 2m | 26CH | 6mm | 100 |
| CT-4030/BULK | Flush cut ends | 3m | 26CH | 6mm | 100 |
| CT-4031/W/BULK | F/M | 3m | 26CH | 6mm | 100 |
| CT-4032/BULK | F/F | 3m | 26CH | 6mm | 100 |
| CT-4034/BULK | F/F/M | 3m | 26CH | 6mm | 100 |
| CT-4034/NS | F/F/M | 3m | 26CH | 6mm | 20 |
| CT-4034/VC/BULK | F/F/M/VC | 3m | 26CH | 6mm | 100 |
| CT-4042/BULK | F/F | 2.7m | 26CH | 6mm | 100 |
| CT-4042/SD/BULK | F/F | 2.7m | 26CH | 6mm | 100 |
| CT-4044/BULK | F/F/M | 2.7m | 26CH | 6mm | 100 |
| CT-4062/BULK | F/F | 6m | 26CH | 6mm | 50 |
| CT-4072/BULK | F/F | 3.5m | 26CH | 6mm | 50 |
| CT-4074/BULK | F/F/M | 3.5m | 26CH | 6mm | 50 |
| CT-4082/BULK | F/F | 1m | 26CH | 6mm | 200 |
| CT-4092/BULK | F/F | 4m | 26CH | 6mm | 50 |
| CT-4094/BULK | F/F/M | 4m | 26CH | 6mm | 50 |
| CT-4312/BULK | F/F | 1.8m | 24CH | 5mm | 100 |
| CT-4313/NS | F/M/VC | 1.5m | 24CH | 5mm | 100 |
| CT-4315/BULK | F/M | 1.8m | 24CH | 5mm | 100 |
| CT-4316/BULK | F/F/VC | 3.5m | 24CH | 5mm | 100 |
| CT-4324/VC/BULK | F/F/M | 2m | 24CH | 5mm | 100 |
| CT-4332/BULK | F/F | 3m | 24CH | 5mm | 100 |

NON-STERILE SUCTION CONNECTING TUBE (CONT)

| Product Code | Type | Length | CH/Tip | Inside Diam. | Qty per Box |
|-----------------|-----------------|--------|--------|--------------|-------------|
| CT-4334/NS | F/F/M | 3m | 24CH | 5mm | 20 |
| CT-4372/BULK | F/F | 3.2m | 24CH | 5mm | 100 |
| CT-4374/BULK | F/F/M | 3.5m | 24CH | 5mm | 50 |
| CT-4622/BULK | F/F | 6m | 30CH | 8mm | 50 |
| CT-4630/BULK | Flush cut ends | 3m | 30CH | 8mm | 50 |
| CT-4632/BULK | F/F | 3m | 30CH | 8mm | 50 |
| CT-4642/BULK | F/F | 4.0m | 30CH | 8mm | 50 |
| CT-4662/BULK | F/F | 6m | 30CH | 8mm | 50 |
| CT-4713/BULK | F/F | 1.8m | 28CH | 7mm | 100 |
| CT-4722/BULK | F/F | 2m | 28CH | 7mm | 100 |
| CT-4724/BULK | F/F/M | 2m | 28CH | 7mm | 100 |
| CT-4732/BULK | F/F | 3m | 28CH | 7mm | 50 |
| CT-4734/BULK | F/F/M | 3m | 28CH | 7mm | 50 |
| CT-4741/VC/BULK | F/F/M | 2.5m | 28CH | 7mm | 50 |
| CT-4743/BULK | F/F | 2.7m | 28CH | 7mm | 50 |
| CT-4772/BULK | F/F | 3.5m | 28CH | 7mm | 50 |
| CT-4781/BULK | F/M | 3.7m | 28CH | 7mm | 50 |
| CT-4783/BULK | F/F | 3.7m | 28CH | 7mm | 50 |
| CT-4792/BULK | F/F | 4m | 28CH | 7mm | 50 |
| CT-4922/BULK | F/F | 2m | 26CH | 6mm | 100 |
| CT-4932/BULK | F/F | 3m | 26CH | 6mm | 100 |
| CT-4934/BULK | F/F/M | 3m | 26CH | 6mm | 100 |
| CT-4936/BULK | F/F | 3m | 26CH | 6mm | 100 |
| CT-4972/BULK | F/F | 3.5m | 26CH | 6mm | 50 |
| CTS-4027/BULK | F/F/M | 2m | 26CH | 6mm | 50 |
| CTS-4032/BULK | F/F | 3m | 26CH | 6mm | 100 |
| CTS-4062/BULK | F/F | 6m | 26CH | 6mm | 50 |
| CTS-4072/BULK | F/F | 3.5m | 26CH | 6mm | 50 |
| CTS-4092/BULK | F/F | 4m | 26CH | 6mm | 50 |
| CTS-4322/NS | F/F | 2m | 24CH | 5mm | 25 |
| CTS-4332/BULK | F/F | 3m | 24CH | 5mm | 100 |
| CTS-4632/BULK | F/F | 3m | 30CH | 8mm | 50 |
| CTS-4640/R/BULK | F/Flush cut end | 3m | 30CH | 8mm | 50 |
| CTS-4642/BULK | F/F | 4m | 30CH | 8mm | 50 |
| CTS30/10/BULK | No Connectors | 10m | 30CH | 30mm | 15 |
| CT-4012/PC | F/F | 1.8m | 26CH | 6mm | 25 |
| CT-4092/PC | F/F | 4.0m | 26CH | 6mm | 15 |
| CT-4074/PC | F/F/M | 3.5m | 26CH | 6mm | 20 |
| CT-4034/PC | F/F/M | 3.0m | 26CH | 6mm | 20 |
| CT-4094/PC | F/F/M | 4.0m | 26CH | 6mm | 15 |
| CT-4743/PC | F/F | 2.7m | 28CH | 7mm | 20 |



Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk



Ivor Shaw Limited t/a Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, UK.
Tel: +44 (0)1332 794880
www.penninehealthcare.co.uk

 @pennine_health

 Pennine Healthcare

 @Pennine_Health