


NON-STERILE PRODUCT RANGE



CONTENTS

About Pennine	3
Pennine Non-Sterile Range	4
<hr/>	
Non-Sterile Link Yankauers	5
Non-Sterile Link Yankauer Sets	6
Non-Sterile Poole Suction Probes	7
Non-Sterile Poole Suction Probe Sets	8
Non-Sterile Micro-Flex Probes	9
Non-Sterile Suction Probes	10
Non-Sterile Suction Probe Sets	11
Non-Sterile Link Bubble Tubing	12



ABOUT PENNINE

Pennine Healthcare is proud to be one of the UK's leading, employee owned, manufacturer and distributor, of single-use sterile and non-sterile Ward & Theatre Consumables and Custom Procedure Packs, as well as one of the leading suppliers of PPE to healthcare settings.

Since our incorporation in 1963, the business has grown year on year, and we market and distribute products in over 50 countries worldwide, employing around 180 people at our purpose-built premises in Derby. The site covers almost 100,000 sq ft, incorporating an area for storage and loading, capable of handling over 4,500 pallets. We are also home to some of the most sophisticated technology in the industry, including three state-of-the-art clean rooms, operating to BS EN ISO 14644-1 standards.

In line with our ISO 13485 accreditation, Pennine has a policy of continual improvement to meet the needs of all our customers, while continual investment in people, training, and technology has enabled us to provide an increasingly wide range of user-friendly, high-specification products. As a healthcare manufacturer, we understand the

importance and requirements of the products we produce, combining this with promoted standards of safety, health, and welfare to provide a healthy and safe working environment for all staff members.

We are committed to minimizing our direct impact on the environment and are continually working towards improving our environmental performance, integrating recognised best practices into all our business operations.

We were proud to become an Employee-Owned Trust in 2021, ensuring the business stayed in the city of Derby, while recognising staff for their hard work and commitment throughout the years.

In 2023, we were delighted to win the prestigious MAKE UK Midlands & East of England Manufacturing Matters Award, in line with our commitment to UK Manufacturing, as well as launching our Commercial Partnership Offering – specialising in Medical Device Packaging, Kit Assembly, and Contract Manufacturing.

PENNINE NON-STERILE PRODUCT RANGE

Our Bulk Non-Sterile range offers superb quality, performance, and solution, to your surgical needs. We strive to achieve the following in our Non-Sterile range:


CHOICE

Whether you are looking for Link Yankauers, Suction Probes, or more, our product range offers a variety of products in a range of lengths, diameters, and connectors, to fulfil yours and your patients' needs and requirements.




SPEED

Our Non-Sterile range provides quicker lead times, reduced packaging, and cheaper costs. Which means you will receive the devices you require, quicker, more cost effective and in a more sustainable manner.






SAFETY

We are fully MDR Certified. Our products are not manufactured using natural rubber latex or phthalates (DEHP), as standard, putting patient safety at the forefront of everything we manufacture.



We are delighted to announce that Pennine Healthcare have been awarded our MDR Certification for Class I Sterile, Class IIa and Procedure Packs.

NON-STERILE LINK YANKAUERS



Single use, rigid suction devices for use with suction connecting tubing for general purpose surgical suctioning.

- Available with or without vacuum control (VC) vent
- Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip


Product Code	Type	Length	CH	Qty per Box
LY-3400/BULK	Short Mini Tip	190mm	15ch	300
LY-3401/BULK	Short Mini Tip VC	190mm	15ch	300
LY-3500/BULK	Mini Fine Tip	255mm	12ch	300
LY-3501/BULK	Mini Fine Tip (VC)	255mm	12ch	300
LY-3600/BULK	Midi Tip	240mm	18ch	300
LY-3600/PR/BULK	Midi Tip (PR)	240mm	18ch	300
LY-3601/BULK	Midi Tip VC	240mm	18ch	300
LY-3700/BULK	Maxi Tip	275mm	22ch	300
LY-3701/BULK	Maxi Tip (VC)	275mm	22ch	300
LY-3800/BULK	Maxi Plus Tip	280mm	28ch	300
LY-3801/BULK	Maxi Plus Tip (VC)	280mm	28ch	300
LY-4400/BULK	Short Mini Super Tip	195mm	15ch	300
LY-4401/BULK	Short Mini Super Tip (VC)	195mm	15ch	300
LY-4420/BULK	Short Mini Super Tip (2 side)	195mm	15ch	300
LY-4500/BULK	Mini Fine Super Tip	255mm	12ch	300
LY-4500/NL/BULK	Mini Fine Super Tip (INT)	255mm	12ch	300
LY-4501/BULK	Mini Fine Super Tip (VC)	255mm	12ch	300
LY-4501/NL/BULK	Mini Fine Super Tip (VC) (INT)	255mm	12ch	300
LY-4600/BULK	Midi Super Tip	240mm	18ch	300
LY-4600/PR/BULK	Midi Super Tip (PR)	240mm	18ch	300
LY-4601/BULK	Midi Super Tip (VC)	240mm	18ch	300
LY-4700/BULK	Maxi Super Tip	275mm	22ch	300
LY-4700/PR/BULK	Maxi Super Tip (PR)	275mm	22ch	300
LY-4701/BULK	Maxi Super Tip (VC)	275mm	22ch	300



Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk

Please note: All dimensions are nominal.
Please note any products featuring a “*” indicates the product is distributed by Pennine.
For Regulatory and Technical information, Technical Data Sheets, are available upon request.

NON-STERILE LINK YANKAUER SETS




Single use, Link Suction Yankauers, pre-connected to Pennine suction connecting tubing in 2m and 3metre lengths.

- Available with or without vacuum control (VC) vent
- Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip
- Tubing available in 2m and 3m lengths

Product Code	Tubing Diam/Length	Qty per box	Yankauer Type
LYS-5520/BULK	6mm / 2m	100	LY-3500 Mini Fine Tip
LYS-5530/BULK	6mm / 3m	50	LY-4400 Short Mini
LYS-5531/BULK	6mm / 3m	50	LY-3501 Mini Fine Tip VC
LYS-5620/BULK	6mm / 2m	100	LY-3600 Midi
LYS-5621/BULK	6mm / 2m	100	LY-3601 Midi VC
LYS-5630/BULK	6mm / 3m	50	LY-3600 Midi
LYS-5631/BULK	6mm / 3m	50	LY-3601 Midi VC
LYS-5720/BULK	6mm / 2m	100	LY-3700 Maxi
LYS-5721/BULK	6mm / 2m	100	LY-3701 Maxi VC
LYS-5730/BULK	6mm / 3m	50	LY-3700 Maxi
LYS-5731/BULK	6mm / 3m	50	LY-3701 Maxi VC
LYS-6620/BULK	7mm / 2.2m	50	LY-3600 Midi
LYS-6830/BULK	7mm / 3m	50	LY-3800 Maxi Plus
LYS-7730/BULK	8mm/ 3m	50	LY-3700 Maxi

NON-STERILE POOLE SUCTION PROBES



Single use, flexible sump suction probe with removable inner suction probe. Specifically for abdominal surgery, reducing trauma to surrounding tissue.

- Available with multiple side eyes to reduce tissue grab and for quick drainage and suctioning
- Removable Inner probe
- Available with straight or angled tip

Product Code	Type	Length mm	CH/Tip	Qty per box
YS-3006/BULK	Flexible straight tip	300	20CH	300
YS-3008/BULK	Flexible (angled tip)	300	20CH	300





Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



We’re here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

NON-STERILE POOLE SUCTION PROBE SETS





Single use, disposable, suction probes pre attached to 2m and 3m connecting tube.

- Probes available 20ch and 30ch
- Pre-attached to 2m and 3m Tubing

Product Code	Tubing Diam/Length	Qty per box	Poole Type
SYS-5920/BULK	6mm x 2m	50	YS-3006 Flexible Straight
SYS-5930/BULK	6mm x 3m	50	YS-3006 Flexible Straight

NON-STERILE MICRO-FLEX PROBES



Single use, suction probes with fingertip vacuum control. Designed for surgery and in applications where fine intricate controlled suction is required.

- Flexible suction tube with fingertip vacuum control
- Available in a range of lengths to suit a wide range of procedures.

Product Code	Length	CH/Tip	Qty per box
YS-3128/BULK	21.3cm	12ch ø 4.00mm	500
YS-3148/14/BULK	22.3cm	14ch ø 4.66mm	500




For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



We’re here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

NON-STERILE SUCTION PROBES




Single use, disposable, suction probes pre attached to 2m and 3m connecting tube.

- Probes available 20ch and 30ch
- Pre-attached to 2m and 3m Tubing

Product Code	Type	Length	CH/Tip	Qty per Box
SP-3800/BULK	Suction probe with male connector	245mm	20CH	300
SP-3900/BULK	Flat Tip Suction Probe with male connector	250mm	20CH	300
SP-4200/BULK	Maxi Suction Probe	280mm	30CH	300
SP-4300/BULK	Flat Tip Maxi Suction Probe	255mm	30CH	300
YS-3000/BULK	Bent Tip Suction Probe	265mm	20CH	300
YS-3001/BULK	Bent Tip Suction Probe (VC)	265mm	20CH	300

*SP-3800/BULK pictured

NON-STERILE SUCTION PROBE SETS



Single use, disposable, suction probes pre attached to 2m and 3m connecting tube.

- Probes available 20ch and 30ch
- Pre-attached to 2m and 3m Tubing

Product Code	CH/Tip	Qty per Box	Tube ID	Probe Type
SPS-5020/W/BULK	20CH	100	6mm	Standard Probe
SPS-5021/BULK	20CH	100	6mm	Standard Probe VC
SPS-6020/BULK	20CH	100	8mm	Maxi Probe

*SP-3800/BULK pictured





Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



We’re here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

NON-STERILE LINK BUBBLE TUBING




Link bubble tubing - available in 30 and 50 metre lengths supplied non-sterile in convenient dispenser bags.

- Tube extruded with a bubble at approx 1m intervals
- Simply cut through the bubble at the required length to create either male or female end which will fit a suction device
- Three diameters available: 5mm, 6mm and 7mm

Product Code	Type	Length	CH/Tip	Inside Diam.	Qty per Box
SL-8005	Plain	30m	24CH	5mm	8 Coils
SL-8055	Plain	50m	24CH	5mm	4 Coils
SL-8356	Plain	50m	26CH	6mm	4 Coils
SL-8007 (International)	Plain	30m	28CH	7mm	7 Coils

TUBING CONNECTOR RIGID




A short length rigid PVC tube that tapers at each end to form male connectors allowing two sections of suction connecting tubing to be joined to increase the length.

- Can be used to extend length of tubing (i.e to connect two tubes together of either the same or different diameters).

Product Code	Qty per Box
DP-501/BULK	500

CAPPED VACUUM CONTROL CONNECTOR



The vacuum control unit is comprised of two PVC male connectors with a hole in between that can be covered by a finger during use to control the vacuum and suction flow rate of a device while in use. The vacuum control unit also has a cap that can be placed to cover the hole to allow for continuous suction without the need of the user to cover the hole with their finger.

- Can be used to connect two devices and allow control of fluid flow. Male ends push fit into suction tubing and/or device. Open cap and use finger pressure to control the vacuum flow.

Product Code	Qty per Box
VC-2/100/BULK	100
VC-2/BULK	1200



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

NON-STERILE SUCTION CONNECTING TUBE



Single Use, Suction connecting tubing to connect suction sources to suction waste collection systems.

- Secure fit to most suction devices and sources
- Available with optional male connector with or without VC

Product Code	Type	Length	CH/Tip	Inside Diam.	Qty per Box
CT-4012/NS	F/F	1.8m	26CH	6mm	25
CT-4014/NS	F/F/M	1.8m	26CH	6mm	25
CT-4020/BULK	Flush cut ends	2m	26CH	6mm	100
CT-4022/BULK	F/F	2m	26CH	6mm	100
CT-4024/BULK	F/F/M	2m	26CH	6mm	100
CT-4030/BULK	Flush cut ends	3m	26CH	6mm	100
CT-4031/W/BULK	F/M	3m	26CH	6mm	100
CT-4032/BULK	F/F	3m	26CH	6mm	100
CT-4034/BULK	F/F/M	3m	26CH	6mm	100
CT-4034/NS	F/F/M	3m	26CH	6mm	20
CT-4034/VC/BULK	F/F/M/VC	3m	26CH	6mm	100
CT-4042/BULK	F/F	2.7m	26CH	6mm	100
CT-4042/SD/BULK	F/F	2.7m	26CH	6mm	100
CT-4044/BULK	F/F/M	2.7m	26CH	6mm	100
CT-4062/BULK	F/F	6m	26CH	6mm	50
CT-4072/BULK	F/F	3.5m	26CH	6mm	50
CT-4074/BULK	F/F/M	3.5m	26CH	6mm	50
CT-4082/BULK	F/F	1m	26CH	6mm	200
CT-4092/BULK	F/F	4m	26CH	6mm	50
CT-4094/BULK	F/F/M	4m	26CH	6mm	50
CT-4312/BULK	F/F	1.8m	24CH	5mm	100
CT-4313/NS	F/M/VC	1.5m	24CH	5mm	100
CT-4315/BULK	F/M	1.8m	24CH	5mm	100
CT-4316/BULK	F/F/VC	3.5m	24CH	5mm	100
CT-4324/VC/BULK	F/F/M	2m	24CH	5mm	100
CT-4332/BULK	F/F	3m	24CH	5mm	100



Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk

Product Code	Type	Length	CH/Tip	Inside Diam.	Qty per Box
CT-4334/NS	F/F/M	3m	24CH	5mm	20
CT-4372/BULK	F/F	3.2m	24CH	5mm	100
CT-4374/BULK	F/F/M	3.5m	24CH	5mm	50
CT-4622/BULK	F/F	6m	30CH	8mm	50
CT-4630/BULK	Flush cut ends	3m	30CH	8mm	50
CT-4632/BULK	F/F	3m	30CH	8mm	50
CT-4642/BULK	F/F	4.0m	30CH	8mm	50
CT-4662/BULK	F/F	6m	30CH	8mm	50
CT-4713/BULK	F/F	1.8m	28CH	7mm	100
CT-4722/BULK	F/F	2m	28CH	7mm	100
CT-4724/BULK	F/F/M	2m	28CH	7mm	100
CT-4732/BULK	F/F	3m	28CH	7mm	50
CT-4734/BULK	F/F/M	3m	28CH	7mm	50
CT-4741/VC/BULK	F/F/M	2.5m	28CH	7mm	50
CT-4743/BULK	F/F	2.7m	28CH	7mm	50
CT-4772/BULK	F/F	3.5m	28CH	7mm	50
CT-4781/BULK	F/M	3.7m	28CH	7mm	50
CT-4783/BULK	F/F	3.7m	28CH	7mm	50
CT-4792/BULK	F/F	4m	28CH	7mm	50
CT-4922/BULK	F/F	2m	26CH	6mm	100
CT-4932/BULK	F/F	3m	26CH	6mm	100
CT-4934/BULK	F/F/M	3m	26CH	6mm	100
CT-4936/BULK	F/F	3m	26CH	6mm	100
CT-4972/BULK	F/F	3.5m	26CH	6mm	50
CTS-4027/BULK	F/F/M	2m	26CH	6mm	50
CTS-4032/BULK	F/F	3m	26CH	6mm	100
CTS-4062/BULK	F/F	6m	26CH	6mm	50
CTS-4072/BULK	F/F	3.5m	26CH	6mm	50
CTS-4092/BULK	F/F	4m	26CH	6mm	50
CTS-4322/NS	F/F	2m	24CH	5mm	25
CTS-4332/BULK	F/F	3m	24CH	5mm	100
CTS-4338/bulk	F/Soft tube end	3m	24CH	5mm	100
CTS-4378/bulk	F/Soft tube end	3.5m	24CH	5mm	100
CTS-4328/bulk	F/Soft tube end	2m	24CH	5mm	100
CTS-4632/BULK	F/F	3m	30CH	8mm	50
CTS-4640/R/BULK	F/Flush cut end	3m	30CH	8mm	50
CTS-4642/BULK	F/F	4m	30CH	8mm	50
CTS30/10/BULK	No Connectors	10m	30CH	30mm	15
CT-4012/PC	F/F	1.8m	26CH	6mm	25
CT-4092/PC	F/F	4.0m	26CH	6mm	15
CT-4074/PC	F/F/M	3.5m	26CH	6mm	20
CT-4034/PC	F/F/M	3.0m	26CH	6mm	20
CT-4094/PC	F/F/M	4.0m	26CH	6mm	15
CT-4743/PC	F/F	2.7m	28CH	7mm	20

Please note: All dimensions are nominal.
Please note any products featuring a “*” indicates the product is distributed by Pennine.
For Regulatory and Technical information, Technical Data Sheets, are available upon request.



Ivor Shaw Limited t/a Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, UK.
Tel: +44 (0)1332 794880
www.penninehealthcare.co.uk



@pennine_health



Pennine Healthcare



@Pennine_Health