



CONTENTS

About Pennine	3
Pennine Surgical Suction Range	4
Suction Connecting Tubing	
Suction Connecting Tubing	5
Soft Suction Connecting Tubing	7
Flexi-Rib Tubing	8
Tubing with Capped Vacuum Control	9
Capped Vacuum Control Connector	9
Tubing Connector Rigid	10
Link Suction Connecting Tubing	11
Link Suction Connecting Tubing (Bubble Tubing)	12
Suction Devices	
Link Yankauers	13
Link Yankauer Sets	15
AirOline Yankauers & Sets	16
Derby Yankauers	17
Selectavac Yankauers	18
Poole Suction Probes	19
Poole Suction Probe Sets	20
Suction Probes	
Suction Probes	21
Suction Probe Sets	22
Micro-Flex Suction Probes	23
Endo-Flex Suction Probes	24
Frazier Probes	25
Easyflow Orthopaedic Suction Set	26
Uterine Aspiration Curette System	
Flexible Uterine Curettes	27
Uterine Aspiration Curettes	28
Wide-Bore Aspiration Tubing (Flexi-Cuff)	29
Accessories	
Camera Drape	30
C-Arm Drape	31



ABOUT PENNINE

Pennine Healthcare is proud to be one of the UK's leading, employee owned, manufacturer and distributor, of single-use sterile and non-sterile Ward & Theatre Consumables and Custom Procedure Packs, as well as one of the leading suppliers of PPE to healthcare settings.

Since our incorporation in 1963, the business has grown year on year, and we market and distribute products in over 50 countries worldwide, employing around 180 people at our purpose-built premises in Derby. The site covers almost 100,000 sq ft, incorporating an area for storage and loading, capable of handling over 4,500 pallets. We are also home to some of the most sophisticated technology in the industry, including three state-of-the-art clean rooms, operating to BS EN ISO 14644-1 standards.

In line with our ISO 13485 accreditation, Pennine has a policy of continual improvement to meet the needs of all our customers, while continual investment in people, training, and technology has enabled us to provide an increasingly wide range of user-friendly, high-specification products.

As a healthcare manufacturer, we understand the importance and requirements of the products we produce, combining this with promoted standards of safety, health, and welfare to provide a healthy and safe working environment for all staff members.

We are committed to minimizing our direct impact on the environment and are continually working towards improving our environmental performance, integrating recognised best practices into all our business operations.

We were proud to become an Employee-Owned Trust in 2021, ensuring the business stayed in the city of Derby, while recognising staff for their hard work and commitment throughout the years.

In 2023, we were delighted to win the prestigious MAKE UK Midlands & East of England Manufacturing Matters Award, in line with our commitment to UK Manufacturing, as well as launching our Commercial Partnership Offering – specialising in Medical Device Packaging, Kit Assembly, and Contract Manufacturing.

PENNINE SURGICAL SUCTION RANGE

Our Surgical Suction Range offers superb quality, performance, and solution, to your surgical suction needs.

We strive to achieve the following in our Surgical Suction Range:


CHOICE

Whether you are looking for tubing, suction devices, probes, or more, our product range offers a variety of products in a range of lengths, diameters and connectors, to fulfil yours and your patients needs and requirements.




SPEED

Our Surgical Suction Range includes colour coding, for easy size identification, whilst certain products are provided in easy to open packaging.







SAFETY

Our products are not manufactured using natural rubber latex or phthalates (DEHP), as standard, putting patient safety at the forefront of everything we manufacture.



In line with our ISO 13485 accreditation, Pennine has a policy of continual improvement in order to meet the needs of the changing international standards for sterile, single use medical devices.

SUCTION CONNECTING TUBING



Single Use, Suction connecting tubing to connect suction sources to suction waste collection systems, suction catheters, yankauers, suction probes and other suction devices.

- Colour coded funnel for diameter identification
- Secure fit to most suction devices and sources
- Available with optional male connector with or without VC - (See codes with suffix M & M/VC)
- Sterile, Single use, Individual peel pouches





Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk

Please note: All dimensions are nominal.
Please note any products featuring a "*" indicates the product is distributed by Pennine.
For Regulatory and Technical information, Technical Data Sheets, are available upon request.

Product Code	Type	Length	CH	Inside Diam. mm	Qty per Box
CT-4311	F/M	1.8m	24CH	5	25
CT-4312	F/F	1.8m	24CH	5	25
CT-4332	F/F	3.0m	24CH	5	20
CT-4335	F/F	3.0m	24CH	5	50
CT-4362	F/F	2.7m	24CH	5	20
CT-4372	F/F	3.5m	24CH	5	20
CT-4314	F/F/M	1.8m	24CH	5	25
CT-4364	F/F/M	2.7m	24CH	5	20
CT-4334	F/F/M	3.0m	24CH	5	20
CT-4374	F/F/M	3.5m	24CH	5	20
CT-4314/VC	F/F/M/VC	1.8m	24CH	5	25
CT-4324/VC	F/F/M/VC	2.0m	24CH	5	25
CT-4334/VC	F/F/M/VC	3.0m	24CH	5	20
CT-4354/VC	F/F/M VC	1.5m	24CH	5	25
CT-4394/VC	F/F/M/VC	2.7m	24CH	5	20
CT-4012	F/F	1.8m	26CH	6	25
CT-4022	F/F	2.0m	26CH	6	25
CT-4042	F/F	2.7m	26CH	6	20
CT-4032	F/F	3.0m	26CH	6	20
CT-4072	F/F	3.5m	26CH	6	20
CT-4092	F/F	4.0m	26CH	6	15
CT-4062	F/F	6.0m	26CH	6	10
CT-4020	F/open end	2.0m	26CH	6	25
CT-4030	F/open end	3.0m	26CH	6	20
CT-4014	F/F/M	1.8m	26CH	6	25
CT-4024	F/F/M	2.0m	26CH	6	25
CT-4044	F/F/M	2.7m	26CH	6	20
CT-4034	F/F/M	3.0m	26CH	6	20
CT-4074	F/F/M	3.5m	26CH	6	20
CT-4064	F/F/M	6.0m	26CH	6	10
CT-4014/VC	F/F/M/VC	1.8m	26CH	6	25
CT-4024/VC	F/F/M/VC	2.0m	26CH	6	25
CT-4034/VC	F/F/M/VC	3.0m	26CH	6	20
CT-4622	F/F	2.0m	30CH	8	10
CT-4632	F/F	3.0m	30CH	8	10
CT-4672	F/F	3.5m	30CH	8	10
CT-4642	F/F	4.0m	30CH	8	10
CT-4624	F/F/M	2.0m	30CH	8	10
CT-4634	F/F/M	3.0m	30CH	8	10
CT-4674	F/F/M	3.5m	30CH	8	10
CT-4634/VC	F/F/M/VC	3.0m	30CH	8	10
* CT-4037	F/F/M	3.0m	26CH	6	50

SOFT SUCTION CONNECTING TUBING



Single Use, Soft Suction Connecting Tubing to connect suction sources to suction waste collection systems, suction catheters, yankauers, suction probes and other suction devices.

- Softer PVC grade tubing for enhanced durability
 - Not to be used in high pressure applications
- Sterile, Single use, Individual peel pouches

Product Code	Type	Length	CH	Inside Diam. mm	Qty per Box
CTS-4338	F/Soft tube end	3m	24CH	5	20
CTS-4378	F/Soft tube end	3.5m	24CH	5	20
CTS-4328	F/Soft tube end	2m	24CH	5	20
CTS-4332	F/F	3.0m	24CH	5	20
CTS-4022	F/F	2.0m	26CH	6	25
CTS-4025	F/M	2.0m	26CH	6	25
CTS-4026	F/M LL	2.0m	26CH	6	50
CTS-4032	F/F	3.0m	26CH	6	20
CTS-4032/WG	F/F	3.0m	26CH	6	20
CTS-4072	F/F	3.5m	26CH	6	20
CTS-4092	F/F	4.0m	26CH	6	15
CTS-4062	F/F	6.0m	26CH	6	10
CTS-4640/R	F/open end	4.0m	30CH	8	10





Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

FLEXI-RIB TUBING





Single use, sterile, Flexi-Rib tubing to connect suction sources to suction waste collection systems, suction catheters, yankauers, suction probes and other suction devices.

- Ribbed connecting tubing for enhanced drapability
- Non-slip outer rib reduces friction
- Resistant to coiling, kinking
- Sterile, Single use, Individual peel pouches
- Available with optional Vacuum Control - (See codes with VC)
- Colour coded funnel for diameter identification

Product Code	Type	Length	CH	Inside Diam. mm	Qty per Box
CT-4922	F/F	2.0m	26CH	6	60
CT-4932	F/F	3.0m	26CH	6	50
CT-4972	F/F	3.5m	26CH	6	50
CT-4924	F/F/M	2.0m	26CH	6	60
CT-4934	F/F/M	3.0m	26CH	6	50
CT-4722	F/F	2.0m	28CH	7	60
CT-4732	F/F	3.0m	28CH	7	50
CT-4772	F/F	3.5m	28CH	7	40
CT-4734	F/F/M	3.0m	28CH	7	50
CT-4774	F/F/M	3.5m	28CH	7	40
CT-4714/VC	F/F/M/VC	1.8m	28CH	7	25
CT-4741/VC	F/F/M/VC	2.5m	28CH	7	20


TUBING WITH CAPPED VACUUM CONTROL



- Sterile double wrapped individual peel pouches

Product Code	Type	Length	CH	Inside Diam. mm	Qty per Box
CT-4015	F/VCU2	1.8m	26CH	6	25
CT-4025	F/VCU2	2.0m	26CH	6	25

CAPPED VACUUM CONTROL CONNECTOR



This Capped Vacuum Control Connector is capped for ease of use to convert a standard suction device into a vacuum controlled suction device. Available sterile and bulk non-sterile.

Product Code	Type	Qty per Box
VCU-2	Sterile	100
VC-2/Bulk	Bulk, non sterile	1200




For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk

TUBING CONNECTOR RIGID





This Tubing Connector Rigid can be used to extend length of tubing, of either the same or different diameters.

- Clear PVC, single use only
- Available sterile and bulk non-sterile
- Designed for use with all types of tubing

Product Code	Type	Qty per Box
DP-501	Sterile	100
DP-501/Bulk	Bulk, non sterile	500

LINK SUCTION CONNECTING TUBING



Single Use, Pre-cut Link bubble connecting tube.

- Pre-cut Link tubing - used for connecting suction sources to suction waste collection systems, suction catheters, yankauer suckers, suction probes and other suction devices.
- Tubing extruded with an integral female/female funnel connections or tapered male/female funnel connections.
- The tube can be modified by cutting with scissors.
- Available in plain PVC
- Sterile, Single use, Individual peel pouches

Product Code	Type	Length	CH	Inside Diam. mm	Qty per Box
CT-4513	Plain F/F	1.8m	24CH	5mm	25
CT-4543	Plain F/F	2.7m	24CH	5mm	20
CT-4511	Plain M/F	1.8m	24CH	5mm	25
CT-4541	Plain M/F	2.7m	24CH	5mm	20
CT-4713	Plain F/F	1.8m	28CH	7mm	25
CT-4783	Plain F/F	3.7m	28CH	7mm	20
CT-4741	Plain M/F	2.7m	28CH	7mm	20
CT-4743	Plain F/F	2.7m	28CH	7mm	20
CT-4711	Plain M/F	1.8m	28CH	7mm	25
AT-4000	F/M	3m	24CH	5mm	25




We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk


LINK CONNECTING TUBING (BUBBLE TUBING)



Our Link Connecting Tubing (Bubble Tubing) is available non-sterile and in 30 and 50 metre lengths, three diameters: 5mm, 6mm and 7mm supplied in convenient dispenser bags. Tube extruded with a bubble at approximately 1m intervals. Simply cut through the bubble at the required length to create either a male or female end to fit a range of suction devices.

Product Code	Type	Length	CH	Inside Diam. mm	Qty per Box
SL-8005	Plain	30m	24CH	5mm	8
SL-8055	Plain	50m	24CH	5mm	4
SL-8056	Plain	50m	26CH	6mm	4
SL-8007	Plain	30m	28CH	7mm	7

LINK YANKAUERS



Our Link Yankauers are rigid, transparent suction devices for use with suction connecting tubing. For general purpose and fine suctioning, as well as suction of large volumes. Also available as a set, pre-attached to 2 or 3 metre tubing.

- Sterile double wrapped in individual peel pouches
- Available with or without VC (vacuum control) vent
- Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip
- Midi Yankauer available with additional Atraumatic Tip to reduce tissue grab

Product Code	Type	Length	CH	Qty per Box
LY-3400	Short Mini Tip	190mm	15CH	100
LY-3401	Short Mini Tip VC	190mm	15CH	100
LY-4400	Short Mini Super Tip	195mm	15CH	100
LY-4401	Short Mini Super Tip VC	195mm	15CH	100
LY-4420	Short Mini Super Tip (2 side)	195mm	15CH	100
LY-3500	Mini Fine Tip	255mm	12CH	50
LY-3501	Mini Fine Tip VC	255mm	12CH	50
LY-4500	Mini Fine Super Tip	255mm	12CH	50
LY-4501	Mini Fine Super Tip VC	255mm	12CH	50
LY-4501/NL	Mini Fine Super Tip VC	255mm	12CH	50
LY-3600	Midi Tip	240mm	18CH	50
LY-3601	Midi Tip VC	240mm	18CH	50
LY-3600/PR	Midi Tip	240mm	18CH	50
LY-4600	Midi Super Tip	240mm	18CH	50
LY-4600/PR	Midi Super Tip	240mm	18CH	50
LY-4601	Midi Super Tip VC	240mm	18CH	50
LY-3700	Maxi Tip	275mm	22CH	50
LY-3700/F	Maxi Tip (with blue connector)	315mm	22CH	50
LY-3700/RF/FF	Maxi Tip (with straight connector)	315mm	22CH	50
LY-3701	Maxi Tip VC	275mm	22CH	50
LY-3702	Maxi Tip (No side eyes)	275mm	22CH	50
LY-4700	Maxi Super Tip	275mm	22CH	50
LY-4700/PR	Maxi Super Tip	275mm	22CH	50
LY-4701	Maxi Super Tip VC	275mm	22CH	50



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk




We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

LINK YANKAUERS (CONTINUED)

Product Code	Type	Length	CH	Qty per Box
LY-3800	Maxi Plus Tip	280mm	28CH	50
LY-3801	Maxi Plus Tip VC	280mm	28CH	50
LY-3803	Maxi Plus Tip VC (No side eyes)	280mm	28CH	50
* AELY3602	Atraumatic tip, soft short	255mm	18CH	50
* AELY3603	Atraumatic tip, soft short, VC	255mm	18CH	50
* AELY3606	Atraumatic tip, short rigid	255mm	18CH	50
* AELY3607	Atraumatic tip, short rigid, VC	255mm	18CH	50
* AELY3608	Atraumatic tip, long soft	263mm	18CH	50
* AELY3609	Atraumatic tip, long soft, VC	263mm	18CH	50
* GC-3600	Midi Super Tip	235mm	18CH	50

LINK YANKAUER SETS



Our Link Yankauer Sets are pre-connected to Pennine Suction Connection Tubing in 2.0 and 3.0 metre lengths. For general purpose surgical suctioning.

- Sterile double wrapped in individual peel pouches
- Available with or without VC (vacuum control) vent

Sets with 2m tubing

Product Code	Tubing Diam. /Length	Qty Per Box	Yankauer Type
LYS-5520	6mm / 2m	10	LY-3500 Mini Fine Tip
LYS-5620	6mm / 2m	10	LY-3600 Midi
LYS-5621	6mm / 2m	10	LY-3601 Midi VC
LYS-5720	6mm / 2m	10	LY-3700 Maxi
LYS-5721	6mm / 2m	10	LY-3701 Maxi VC
LYS-6620	7mm / 2.2m	10	LY-3600 Midi

Sets with 3m tubing

Product Code	Tubing Diam. /Length	Qty Per Box	Yankauer Type
LYS-5430	6mm / 3m	10	LY-3400 Short Mini
LYS-5530	6mm / 3m	10	LY-3500 Mini Fine Tip
LYS-5531	6mm / 3m	10	LY-3501 Mini Fine Tip VC
LYS-5630	6mm / 3m	10	LY-3600 Midi
LYS-5631	6mm / 3m	10	LY-3601 Midi VC
LYS-5730	6mm / 3m	10	LY-3700 Maxi
LYS-5731	6mm / 3m	10	LY-3701 Maxi VC
LYS-6830	7mm / 3m	5	LY-3800 Maxi Plus
LYS-7730	8mm / 3m	5	LY-3700 Maxi




Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk

STERILE AIROLINE YANKAUERS & SETS





Single use, rigid suction devices for use with suction connecting tubing for general purpose surgical suctioning. AirOline set, yankauer pre-connected to Pennine suction connecting tubing in 1m and 2m lengths.

- Available as a Maxi Plus Tip with a smooth bend
- The Pennine AirOline Yankauer devices are designed with a smooth bend to help with insertion into the patient airway.
- This device is appropriate for the management of difficult airways.
- Tubing available in 1m and 2m lengths
- This device can aid in the management of a patient's airway during an emergency event.
- Sterile, double wrapped individual peel pouches

Product Code	Product Description	Tubing Diameter / Length	Yankauer Type	Yankauer Length / CH	Qty per box
LY-3800/DB	Pennine AirOline Yankauer Maxi Plus Tip 28fg 265mm	N/A	Maxi Plus Tip	265 mm / 28 CH	50
LYS-5810/WW/DB	Pennine AirOline Yankauer Set, Tube 6mm x 1m, Maxi Plus Tip 28fg 265mm	6mm / 1m	Maxi Plus Tip	265 mm / 28 CH	15
LYS-5820/WW/DB	Pennine AirOline Yankauer Set, Tube 6mm x 2m, Maxi Plus Tip 28fg 265mm	6mm / 2m	Maxi Plus Tip	265 mm / 28 CH	10

DERBY YANKAUERS



Our Derby Yankauers are clear, acrylic suction devices, similar to traditional yankauer suction probes. With or without sump tip to reduce potential tissue suction during use.

- Sterile double wrapped in individual peel pouches
- Available with or without VC (Vacuum control) vent
- Available with or without sump tip to reduce the risk of tissue damage
- Ergonomically shaped handle for a comfortable grip during use

Product Code	Type	Length	Diam.	CH	Qty per Box
* AEYS3004	with sump/rose tip	270mm	4.2mm	24CH	50
* AEYS3005	with sump/rose tip, VC	270mm	4.2mm	24CH	50
* AEYS3014	with straight/flat tip	270mm	4.2mm	24CH	50
* AEYS3015	with straight/flat tip, VC	270mm	4.2mm	24CH	50
YS-3000	Midi Yankauer with straight tip	265mm	6.66mm	20CH	50
YS-3001	Midi Yankauer with VC	265mm	6.66mm	20CH	50





Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk

SELECTAVAC YANKAUERS



Our Selectavac Yankauers are surgical suction probe with switch on and off valve, available with bulb tip and straight tip.

- Switch system allows the operator to manually open and close suction
- Supplied sterile, double wrapped in individual peel pouches
- Easy connection to Pennine range of tubing

Product Code	Type	Length	Diam.	CH/Tip	Qty per Box
* AEYS3009	Bulb Tip, with On/Off Slide Switch	265mm	4.2mm	24CH	50
* AEYS3010	Straight Tip, with On/Off Slide Switch	265mm	4.2mm	24CH	50

POOLE SUCTION PROBES






AEYS3007



YS-3006



YS-3008



Single use, flexible sump suction probe with removable inner suction probe. Specifically for abdominal surgery.

- Available with angled tip
- Sterile double wrapped in individual peel pouches
- Removable inner probe

Product Code	Type	Length	Inside Diam. mm	CH/Tip	Qty per Box
YS-3006	Flexible straight tip	300mm	6.66	20CH	50
YS-3008	Flexible angled tip	300mm	6.66	20CH	50
* AEYS3007	Rigid (white)	264mm	6.5	15CH	50




Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk



Want more information?
Get in touch on +44 (0) 1332 794880 or email customer.services@penninehealthcare.co.uk

POOLE SUCTION SETS




Single use, flexible ‘sump’ suction probe pre-attached to 2m connecting tube. For deep abdominal suctioning or where a ‘sump’ suction action is required.

- Sterile double wrapped in individual peel pouches

Set with 2m tubing	CH/Tip	Qty per Box	Tube Diam. mm	Poole Type
SYS-5920	20CH	15	6	YS-3006
SYS-5923	20CH	15	6	YS-3008
SYS-5930	20CH	10	6	YS-3006

SUCTION PROBES



Our Suction Probes are manufactured from rigid transparent PVC, for general suction purposes in a wide variety of procedures. Easily attached to Pennine suction tubing.

- Sterile double wrapped in individual peel pouches
- Standard probes 20CH

Product Code	Type	Length	CH/Tip	Qty per Box
MP-3802	Short Suction Probe (no side eyes)	120mm	20CH	100
SP-3800	Suction probe with male connector	245mm	20CH	50
SP-3900	Flat tip probe with male connector	250mm	20CH	50
SP-4200	Maxi Suction Probe	280mm	30CH	50
SP-4300	Flat Tip Maxi Suction Probe	255mm	30CH	50





For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk

SUCTION PROBE SETS



Our Suction Probe Sets are available pre-attached to 2m or 3m tubing. For general suction purposes in a wide variety of procedures.

- Sterile double wrapped in individual peel pouches
- Standard probes 20CH

Sets with 2m tubing

Product Code	CH/Tip	Qty Per Box	Tube Diameter	Probe Type
SPS-5020/G	20CH	25	6mm	Standard Probe
SPS-5020/W	20CH	25	6mm	Standard Probe
SPS-5021	20CH	25	6mm	Standard Probe VC
SPS-6020	30CH	10	8mm	Maxi Probe



Sets with 2.1m tubing

Product Code	CH/Tip	Qty Per Box	Tube Diameter	Probe Type
SPS-5022	20CH	50	6mm	Standard Probe

Sets with 3m tubing

Product Code	CH/Tip	Qty Per Box	Tube Diameter	Probe Type
SPS-5030W	20CH	20	6mm	Standard Probe

MICRO-FLEX SUCTION PROBES



Our Micro-Flex Suction Probes are made from PVC and have a fingertip vacuum control. Designed for fine suctioning and in applications where controlled suction is required.

- Sterile double wrapped in individual peel pouches

Product Code	Length	CH/Tip	Qty per Box
YS-3068	16.8cm	6ch ø 2.00mm	50
YS-3088/Z	16.8cm	8ch ø 2.66mm	50
YS-3088	22.3cm	8ch ø 2.66mm	50
YS-3108	22.3cm	10ch ø 3.33mm	50
YS-3128	21.3cm	12ch ø 4.00mm	50
YS-3148/14	22.3cm	14ch ø 4.66mm	50





For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk

ENDO-FLEX SUCTION PROBES






These Endo-Flex Suction Probes are used for general suction where longer probe length is required. All 12CH tubes with varying probe length. Fingertip vacuum control with tapered male connector end of prestrut to fit Pennine connection tubing.

- Sterile double wrapped in individual peel pouches

Product Code	Length	CH/Tip	Qty per Box
YS-3128/20	29.8cm	12ch ø 4.00mm	10
YS-3128/26	35.8cm	12ch ø 4.00mm	10
YS-3128/40	49.8cm	12ch ø 4.00mm	10
YS-3128/60	69.8cm	12ch ø 4.00mm	10

FRAZIER PROBES



Our Frazier Probes are sterile single-use probes. With 'VC' vacuum control hole for precise control over the amount of suctioning. Used in open surgical procedures where delicate suctioning is required.

- Malleable shaft can be bent to desired position
- Sterile double wrapped in individual peel pouches
- Longer probe
- Longer, non-slip handle

Product Code	Length	CH/Tip	Qty per Box
* AEMF1108	230mm	8CH	40
* AEMF1110	230mm	10CH	40
* AEMF1112	230mm	12CH	40




We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk

EASYFLOW ORTHOPAEDIC SUCTION SET




This Easyflow Orthopaedic Suction Set is a sterile filtered suction system for orthopaedic surgery.

- Suitable for removal of large bone fragments during orthopaedic surgery
 - Eliminates blockages
 - Rapid filter exchange
- Supplied sterile
 - Ergonomic handle design ensures comfort and control
 - Reversible suction tip provides extra length for intra-femoral suction

Product Code	Use	Contents	Qty per Box
* AEEFS1111	Orthopaedic Surgery	- 7mm/3m Flexi-Rib Tubing - 2 curved probes - 1 straight probe - Spare filter	10

FLEXIBLE UTERINE CURETTES



Uterine Suction Catheter connected to uterine suction tube is used for curettage and evacuation of the endometrial contents during induced termination of pregnancy or after miscarriage or as a curative measure in cases of dysmenorrhea.

- Printed markings for sizing and correct tip orientation in Flexible Currettes
 - Sterile, double wrapped, individual peel pouches
- Leak proof seals

Product Code	Type	Diam. mm	Qty per Box
* AEAC3004	Flexible	4mm	50
* AEAC3005	Flexible	5mm	50
* AEAC3006	Flexible	6mm	50
* AEAC3007	Flexible	7mm	50
* AEAC3008	Flexible	8mm	50
* AEAC3009	Flexible	9mm	50
* AEAC3010	Flexible	10mm	50
* AEAC3011	Flexible	11mm	50
* AEAC3012	Flexible	12mm	50




For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk

UTERINE ASPIRATION CURETTES




Uterine Suction Catheter connected to uterine suction tube is used for curettage and evacuation of the endometrial contents during induced termination of pregnancy or after miscarriage or as a curative measure in cases of dysmenorrhea.


- Curved and straight design for options in Rigid Curettes
- Moulded markings for sizing and correct tip orientation in Rigid Curettes
- Sterile, double wrapped, individual peel pouches
- Leak proof seals

Product Code	Type	Diam. mm	Qty per Box
* AEKC1506	Straight	6mm	50
* AEKC1507	Straight	7mm	50
* AEKC1508	Straight	8mm	50
* AEKC1509	Straight	9mm	50
* AEKC1510	Straight	10mm	50
* AEKC1511	Straight	11mm	50
* AEKC1512	Straight	12mm	50
* AEKC1006	Curved	6mm	50
* AEKC1007	Curved	7mm	50
* AEKC1008	Curved	8mm	50
* AEKC1009	Curved	9mm	50
* AEKC1010	Curved	10mm	50
* AEKC1011	Curved	11mm	50
* AEKC1012	Curved	12mm	50
* AEKC1514	Straight, with connector	14mm	30
* AEKC1014	Curved, with connector	14mm	30

WIDE-BORE ASPIRATION TUBING (FLEXI-CUFF)



*AEKHS3000M pictured



This Wide-Bore Aspiration Tubing (Flexi-Cuff) offers a convenient and adaptable system for uterine aspiration. It is frequently connected directly or through an external tube or catheter to an irrigator and/or aspirator for irrigation, infusion, or withdrawal of substances.

- Sterile double wrapped individual peel pouches
- Loose male connector in each pack for connection to suction machines where required
- Leak proof seals
- Tubing made of PVC
- Easily connects to curettes
- Flexi-Cuff provides extra flexibility at a crucial point offering ease of use and improved technique

Product Code	Type	Tubing Length	Qty per Box
* AEKHS3000	Wide Bore Tubing + Rotary Handle, loose male connector	2.0m	10
* AEKHS3000M	Wide Bore Tubing + Rotary Handle, attached male connector	2.0m	10
* AEKHS3030	Wide Bore Tubing + Rotary Handle, loose male connector	3.0m	10




For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk



For more information on our product range please contact us on +44 (0) 1332 794880 or alternatively send an email to customer.services@penninehealthcare.co.uk

CAMERA COVER




The disposable camera sleeve is a sterile, single-use disposable drape made of tough, thin polyethylene. The drape is the quickest and simplest way to create a sterile field for use with video/ camera systems during operative laparoscopic procedures. The tape used to secure the camera drapes is elastic and allows a secure and tight closure.

- Made of high grade special anti-blocking PE film
- One strip of adhesive tape
- Film stored in a ring.
- Supplied sterile

Product Code	Description	Size	Qty per Box
* E6222	Camera Cover, Open Tip, Ring Applicator	249cm x 15cm	64

C-ARM DRAPE



*E6018 product drawing pictured

This 3-piece C-Arm drapes is designed to cover x-ray image intensifier equipment and protect the sterile field during surgical procedures.

- Two bonnet covers for the x-ray tank and image intensifier head
- Integral clips for secure fit

Product Code	Description	Size	Qty per Box
* E6016	C-Arm Kit V	76cm x 91cm, 224cm	20
* E6017	C-Arm Kit VI	79cm x 91cm, 224cm	20
* E6018	C-Arm Kit VII	76cm x 91cm, 218cm	20



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk



We're here to help, please get in touch via telephone on +44 (0) 1332 794880 or email us at customer.services@penninehealthcare.co.uk



Ivor Shaw Limited t/a Pennine Healthcare
300 City Gate, London Road, Derby, DE24 8WY, UK.
Tel: +44 (0)1332 794880
www.penninehealthcare.co.uk

 @Pennine_Health  Pennine Healthcare  @Pennine_Health

BRO-7_Rev_13 10/2025