



It promotes a standardised approach to practice, aligns with Difficult Airway Society (DAS) guidelines, and helps improve patient safety when it matters most.

Although eFONA is rarely performed, a 'can't intubate, can't oxygenate' (CICO) scenario is a high-risk, time-critical emergency where preparation is essential. Having a dedicated, standardised procedure pack ensures immediate access to the required equipment, reduces variation in clinical practice, and supports fast, effective intervention when every second counts.



### What Components Does Your eFONA (Emergency Front of Neck Airway) Procedure Pack Include?

- ✓ Scalpel, size 10, safety
- ✓ Swabs, 10cm x 10cm
- ✓ Angled Bougie, 15ch x 600m
- ✓ 6.0mm, cuffed Endotracheal tube
- ✓ Lubricating Jelly Sachet
- ✓ Syringe, 20ml, luer slip
- ✓ Cotton tie
- ✓ Tray

### Key Benefits Of Our eFONA (Emergency Front of Neck Airway) Procedure Pack:

- ✓ Rapid access to essential equipment in emergency situations
- ✓ Supports compliance with Difficult Airway Society (DAS) guidelines
- ✓ Promotes a standardised clinical approach
- ✓ Enhances patient safety
- ✓ Reduces cognitive load in high-pressure scenarios

# WANT TO FIND OUT MORE?

To learn more about the latest addition to our packs portfolio, get in touch now!



[www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)



[customer.services@penninehealthcare.co.uk](mailto:customer.services@penninehealthcare.co.uk)



+44 (0)1332 794880